



Vital Sense O₂ (10L) Oxygen Concentrator

SKU: 83010

The VitalSense O₂ is a high performance medical device that extracts oxygen from atmospheric air to provide reliable supplemental oxygen therapy. Designed to supply patients with a steady oxygen flow of up to 10L/minute, it ensures consistent therapy with an oxygen concentration of 93% $\pm$ 3%.



Flowmeter



LED Display



Humidifier Bottle



Features

- Delivers 93% $\pm$ 3% oxygen at a range of 1-10L/min for reliable therapeutic use
- Operates at 50dB(A), ensuring minimal noise for home or facility environments
- High pressure/High temperature Safety shutdown
- Powered by a standard AC120V outlet with an 800VA power consumption for high-capacity performance
- Precise adjustable oxygen flow control from 1 to 10L/min
- Maintains a safe and consistent output pressure of 20kPa–60kPa
- Audible & Visual O₂ & System Failure Alarms
- Produces high-purity medical grade oxygen
- Utilizes Molecular Sieve (Artificial Zeolite) to efficiently separate oxygen from ambient air
- Compatible with nasal cannulas and humidifier bottles

Specs

Item #	Dimensions	Weight	Warranty
83010	13.78" x 9.84" x 26.38"	56 lb	2 Years or 5,000hrs

