

# Easy Rise *Max* Patient Lift

Experience safe, smooth patient transfers without limits with our full electric Easy Rise Max Patient Lift. It boasts a robust 1000 lbs weight capacity ensuring you meet every patients need, no matter their size. Like our other lifts it too features a convenient battery swap capability to ensure uninterrupted operation.

**1000 LBS  
Weight  
Capacity**

**Care Without Limits**  
*Embrace the potential of  
limitless care with the Easy  
Rise Max.*

## Features

- 6-point spreader bar
- 1000 lbs weight capacity
- Electric powerbase
- Facilitates safe and smooth patient transfers
- Standalone battery option for added convenience



- Intuitive controls for easy and precise lifting and lowering
- Quick and hassle-free attachment of slings
- 360-degree swivel for seamless patient positioning
- Seamless battery swap capability for uninterrupted use



# Spec Table

## General Specifications

Title	Value
Model #	71000
Max. Hanger Bar Height	76.5"
Min. Hanger Bar Height	36.5"
Base Width (Open)	56.4"
Base Width (Closed)	32.6"
Base Height (Clearance)	4.5"
Base Length	60.5"
Caster Size (Front/Rear)	3" / 5"
Weight Capacity	1000 LBS

## Electrical Specifications

Title	Value
Battery	24V DC
Charger Input	100-240V AC
Accessories (Optional)	Standalone Battery Charger
Battery Status	LCD Display

