



Expandable Bariatric Bed



The Infinity bed is a fully electric bed designed with safety in mind. With integrated width and length expandability, a 650 lb. capacity, low position, and many more advanced features, the Infinity bed will become every facility's preferred bed by residents and staff alike.

SPECS

Model	54000
Bed Width Range	Integrated 36" - 42" - 48"
Bed Length Range	Integrated 80" - 84" - 88"
Bed Height Range	7" - 30"
Weight Capacity	650 lbs.
Warranty	15 yrs. Welds/ 5 yrs. Electric/ 1 yr Components



INTEGRATED
EXPANSION



COMFORT CHAIR
POSITION



MOBILITY AT
ANY HEIGHT

FEATURES

- Integrated width/length expansion
- Comfort chair position
- Auto contour
- Battery backup
- Mobility at any height
- Central locking system
- Optional USB port
- Quick release head & footboards
- Underbed lighting
- Bumper guard
- Safety stop
- Heel lift extension
- Autoglide



Expandable Bariatric Bed

✉ sales@emeraldsupply.net | 🌐 emeraldsupply.net



Bumper guard



Integrated length expansion



One touch central lock



Integrated width expansion




Auto Contour



Comfort chair position



Under bed light



RAILS



Q Assist Bar



Pivot Assist Bar



Soft Style Side Rails

BARS & RAILS	ITEM #
Soft Style Side Rails	52210
Q Assist Bar	52211
Pivot Assist Bar	52212

HEAD & FOOTBOARDS



Storm Gray



Solar Oak



Mahogany



Amber Cherry



Expandable Boards

BOARD COLORS	ITEM #
Solar Oak	54020
Mahogany	54021
Amber Cherry	54022
Storm Gray	54023