



Easy Rise 700 Patient Lift

**700 LBS
Weight
Capacity**

Experience safe, smooth patient transfers with our fully electric Easy Rise Patient Lift. It features a robust 700 lbs weight capacity helping you accommodate a wider variety of patients, and its convenient battery swap capability ensures uninterrupted operation.



Features

- Electric Powerbase
- 6-point spreader bar
- 700 lbs weight capacity
- Facilitates safe and smooth patient transfers
- Standalone battery option available for added convenience
- Intuitive controls for easy and precise lifting and lowering
- Quick and hassle-free attachment of slings
- 360-degree swivel for seamless patient positioning
- Seamless battery swap capability for uninterrupted use



Spec Table

General Specifications

Title	Value
Model #	71700
Max. Hanger Bar Height	74.7"
Min. Hanger Bar Height	15.6"
Base Width (Open)	49.1"
Base Width (Closed)	28"
Base Height (Clearance)	4.1"
Base Length	50.7"
Caster Size (Front/Rear)	3" / 4"
Weight Capacity	700 LBS

Electrical Specifications

Title	Value
Battery	24V DC
Charger Input	100-240V AC
Accessories (Optional)	Standalone Battery Charger
Battery Status	LCD Display

