

# Gel Infused Pressure Reducing Mattress

The gel infused foam mattress reduces the likelihood of pressure injuries, and creates a cooling effect channeling heat away from the body, regulating temperature, and promoting airflow. The gel foam redistributes pressure across the surface of the mattress, which leads to the reduction of cradling at higher levels of compression. Thoughtfully designed for both comfort and function, this mattress provides a stable and supportive surface while helping to reduce pressure and improve overall sleep quality.



## Features

- Gel infused foam
- Cooler sleep surface
- 4 Layer Foam Construction
- Optimal Pressure redistribution over large surface area
- High-density foam perimeter provides support and safe patient transfer
- Fluid resistant, easy to clean, vapor permeable stretch mattress cover, ensuring maximum user comfort

## Specs

Item #	65000
Size	42" x 80" x 6"
Weight Capacity	650 lbs.
Warranty	15Yr Core / 2Yr Cover

