

OASIS EXPANDABLE

Hi-Lo Electric Bed



The Oasis Expandable Bed is designed to deliver superior comfort, durability, and safety, with enhanced features for those requiring additional support. Unlike the standard Oasis, this model offers a higher weight capacity and an integrated width expandability of up to 48 inches, accommodating a broader range of patient needs. Its adjustable frame and ultra-low positioning provides increased safety and user flexibility.

SPECS

Bed Length Range	Integrated 80" - 84" - 88"
Bed Width Range	Integrated 36" - 42" - 48"
Bed Height Range	7" - 30"
Weight Capacity	550 lbs
Model	52300

FEATURES


- Length & width expansion
- Expandable head & footboards
- Mobility at any height
- Central locking system
- 8 function hand control
- Heel lift extension
- Bumper guard
- Auto contour
- 2 caster guide locks
- Optional USB port
- Autoglide



One touch central lock



Integrated width expansion



Integrated length expansion



8 function hand control

RAILS

BARS & RAILS	ITEM #
Soft Style Side Rails	52210
Q Assist Bar	52211
Pivot Assist Bar	52212



Q Assist Bar



Pivot Assist Bar



Soft Style Side Rails

HEAD & FOOTBOARDS

BOARD COLORS	ITEM #
Solar Oak	54020
Mahogany	54021
Amber Cherry	54022



Solar Oak



Mahogany



Amber Cherry



Expandable Boards