




## Comfort Memory+ Viscoelastic Infused Mattress

Engineered for clinical grade support, the Comfort Memory Plus Foam Mattress combines advanced viscoelastic memory foam with a multi-layer design to deliver exceptional comfort and durability. It gently contours to the users body, relieving pressure points and promoting proper spinal alignment for a healthier sleep posture. The balanced foam layers provide responsive support while minimizing pressure on any given area.


Model	61091
Dimensions	42" x 80" x 6"
Weight Capacity	450 lbs.
Foam Layers	4
Fire Rating	CFR 1632/CFR 1633
Non-Prorated Warranty	15 Year Core / 3 Year Cover

### Key Features

- Viscoelastic memory foam top layer
- High resilience horizontal die cut foam
- Pressure offloading sloped heel section
- Reinforced Foam Perimeter



## NEW Cover Improvements



Fluid resistant stretch cover remains highly breathable, but the seams are now welded; reduced sew holes means higher water resistance



Easy to use and intuitive tie-down straps quickly and easily secure mattress to virtually any bed frame to decrease mattress migration and reduce falls



All mattresses now include upgraded handles for easier moving and positioning



Upgraded fabric changed to weldable PU-coated ticking for softer, more supple, longer-lasting cover fabric

